



**SB ARTS AND KCP SCIENCE COLLEGE
VIJAYAPURA**

SEMINAR

PLANKS CONSTANT

BY

SAMEER WALIKAR

04

S. B. Walikar

SADDEEQ WALIKAR

05

S. B. Walikar

SOHAIL KHATIB

03

S. B. Khatib

VIKRAM

02

V. K. Ram

SHAMEER JATGAR


18

S. B. Jatgar


AFRIDI NABIWALE

20

A. Fridi




“When you change the way you look at things, the things you look at change.”



TABULAR COLUMN AND NATURE OF GRAPH


| COLOUR OF LED | WAVELENGTH λ in nm | $\frac{1}{\lambda} \cdot 10^{-3}$ | KNEE VOLTAGE V | λV | $h = \frac{e\lambda V}{c}$ In js |
|---------------|----------------------------|-----------------------------------|----------------|-------------|-------------------------------------|
| red | 660 | | | | |
| yellow | 580 | | | | |
| green | 540 | | | | |
| white | 550 | | | | |
| orange | 590 | | | | |

Knee voltage is the voltage at which the bulb start glowing .



Aim : Determination of Planck's

- **Apparatus :** power supply, a one way key, a rheostat, a digital milliammeter, a digital voltmeter, a 1K resistor and different known wavelength LED's (Light-Emitting Diodes).




- **Formula** Planck relation:

$$E = hv = \frac{hc}{\lambda}$$

where:

E = energy
h = Plank constant
v = frequency
c = speed of light
 λ = wavelength

CONSTRUCTION AND PROCEDURE



CALCULATIONS

$$h = \frac{e\lambda v}{c}$$

$$c = 3 \times 10^8 \text{ m/s}$$

$$e = 1.60217662 \times 10^{-19} \text{ coulombs}$$

result:

hence we determine the value of planks constant of different color's of different wavelength by taking average of all the results.



BLDEA's Association

SB ARTS & KCP SCIENCE COLLEGE
VIJAYAPUR -586103

Subject : HALL EFFECT

Presented By :

| Name : | No. |
|----------------|------------|
| Ayesha | 06 |
| Nikita | 10 |
| Megha | 13 |
| Ashwini | 19 |
| Sonali | 16 |

INTRODUCTION:

Hall effect

Definition :

If a piece of conductor (metal or semiconductor) carrying current is placed in a transverse magnetic field, an electric field is produced inside the conductor in a direction normal to the both the current and the magnetic field. This phenomenon is known as the Hall Effect. And the generated voltage is called as the Hall Voltage.

Aim:

Use of Hall effect to setup determine the Hall coefficient for a given semiconductor specimen

1. Hall coefficient
2. Hall mobility
3. Carrier concentration

Apparatus:

Hall probe, connecting wires, gauss meter, digital voltmeter, 2 solenoids, etc.

FORMULA:

Hall coefficient $RH = \frac{V_H}{HI} (\ln \Omega M/T) = \frac{t^*}{H}$ Slope

Carries concentration $n = 3 \frac{\pi}{8eRH} (\text{in } m^{-3})$

Hall mobility $\mu = \frac{8}{3\pi} \sigma RH = \frac{8}{3\pi} RH * \frac{1}{\rho} (\text{in } m^{-2} s^{-1})$

where $\rho = \frac{bt^* V_H}{I}$

Relaxation time, $\tau = \frac{1}{\rho} (\text{in } \rho)$

where V_H -Hall voltage (in volts)

I-Magnetic field (in A)

μ -Magnetic field (in tesla)

E-charge of electron

m-mass of electron

t-thickness of specimen (in m)

l-length of specimen (in m)

b-breadth of specimen (in m)

Circuit Diagram:






Fig. 10.42 Set-up for the measurement of Hall voltage



Nature of graph:



Theory:

Working:

A rectangular slab of the given material having a thickness t and width b is taken and a current I_x ampere is allowed to pass through this sample by connecting it to the battery B_1 . The slab is then placed between two pole pieces of an electromagnet such that if the current direction coincides with the x -axis, the magnetic flux density B coincides with x -axis.

The Hall voltage V_H is then measured by placing two probes at the centres of bottom and top faces of the sample. If the magnetic flux density is B weber/m² and the Hall voltage is V_H volts, then R_H is obtained.

Procedure

- Connect constant current source to the solenoid
- Probes are connected to the gauss meter and placed at the middle of the two solenoids.
- Switch on the gauss meter and constant current source.
- Vary the current through the solenoid from 1: mp to 5 Amp with the interval of 0.5 A, and note the corresponding gauss meter readings.
- Switch off the gauss meter and constant current source and turn the knob of constant current source towards min current
- Fix the hall probe on a wooden stand. Connect green wire to constant current generator and connect red wires to millivoltmeter in the hall effect apparatus
- Replace the probe with hall probe and place the sample material at the middle of the two solenoids
- Switch on the constant current source and CCG
- Carefully increase the current I from CCG and measure the corresponding hall voltage V_H . Repeat this step for different magnetic field B
- Thickness t of the sample is measured using screw gauge.
- Hence calculate the hall coefficient and carrier concentration "n."

Observation:

1. Sample used-Germanium
2. Thickness of sample $t=0.5\text{mm}$
3. Length of sample $l=7\text{mm}$
4. Breadth of sample $b=7\text{mm}$
5. Charge of electron $e=1.6*10^{-19}$
6. Mass of electron $m=9.1*10^{-31}\text{kg}$



BLDEA' s Association

SB ARTS & KCP SCIENCE PU COLLEGE VIJAYAPUR -586103

Subject : Inverse Square Law

Presented By :

Name :

Ayesha

Nikita

Megha

Ashwini

Sonali

Roll No.

06 *Anubandary*


10 *AN*

13 *Megha*

19 *Ashwini*

16 *Sonali*

Circuit Diagram



PROCEDURE:

- Connect the gm to counter scale and power supply.
- Take 3 Blank readings using counter & mean of them gives the background reading B.
- Place the source on the holder at a fixed distance D from the tube take two readings and calculate the mean of them.
- Subtract b from mean gives corrected count per 30 sec & count per second gives intensity I of radiation reaching the tube with distance D between them calculate $1/\text{square root of I}$
- Now change the distance D between source & tube continue step 3 for about 8-10 10 values of D
- Plot graph of $1/\text{square root I V/S D}$